



UltraSal-2

Oral Fluid Collection Device

Oral Fluid is rapidly becoming the sample of choice for a wide range of applications. It's readily available, easy to collect and the collection can be monitored.

UltraSal-2 is ideal for collecting oral fluid samples to test for the presence of abused drugs, hormones, pharmaceuticals, pathogens, genetic markers and various research applications.

The UltraSal-2 offers the following advantages:

- Neat Saliva Collection
- Graduated Volume Indicator
- Split-Sample Collection
- Ideal for Research Applications



Neogen also offers the NeoSal™ Oral Fluid Collection System which features an extraction buffer and collection pad with a built-in sample volume adequacy indicator. The NeoSal is intended for the collection and transportation of oral fluid specimens and is for Forensic Use Only.



Ordering Information	
Product No.	Description
US2-50	UltraSal-2 pack of 50 devices Specifications: Collection tubes: 17 mm x 100 mm; caps: 17 mm Simple assembly required. Snap collection tubes and mouthpiece to the base.
128101-25	NeoSal pack of 25 devices



800-477-8201 (USA/Canada) • 859-254-1221
inform@neogen.com • toxicology.neogen.com

